

Déclaration d'intégration IHE			
Vendeur	Product Name	Version	Date
EOS imaging	EOS_34	1.7.5	12/03/15

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Consistent Time	Time Client	None
Radiation Exposure Monitoring	Acquisition Modality	None
Scheduled Workflow.b	Acquisition Modality	Broad Worklist Query
Scheduled Workflow.b	Acquisition Modality	None

Internet address for vendors IHE information :	
Links to Standards Conformance Statements for the Implementation	
HL7	
DICOM	
Link to general information on IHE : www.ihe.net	